


## Medical/Bio Research Topics I : Week 04 (25 March 2025)

# **Structural MRI (2): Data Processing Methods**

**구조 자기공명영상 (2):  
데이터 처리 방법**



# Brain Mapping with Structural MRI (sMRI)

- T1/T2-weighted sMRI




# Automated Abnormality Detection

- Leverages deep learning algorithms to detect pathological changes in brain tissue
  - Identifies subtle abnormalities that may be overlooked in visual assessment
- Segmentation: Automated delineation of lesion boundaries with voxel-level precision
- Grading: Standardized classification of abnormality severity using predefined scales



[Tomita et al., 2020]

## 3D Convolutional Neural Network for Stroke Lesion Segmentation




[Wang et al., 2024; <https://github.com/mileswyn/SAMIHS>]

## Segment Anything Model (SAM)-based Fine-Tuning for Stroke Lesion Segmentation


# Preprocessing

- Numerous steps to clean and standardize sMRI data before brain morphometry
  - Correction for bias fields (intensity non-uniformities)
    - From a broader range of sources, including imperfections in the MRI scanner's main magnetic field, inhomogeneities in the radiofrequency coil performance, and magnetic susceptibility-induced field inhomogeneities
    - Often characterized by a smooth variation in image brightness

- Segmentation
  - Classifies an image into the non-brain and brain and, furthermore, the brain into different tissues usually including grey matter, white matter, and cerebrospinal fluid
- Normalization
  - Transforms an image from a native space to the standard space usually in the Montreal Neurological Institute (MNI) coordinate system



Intensity non-uniformity



**Correction for Intensity Non-uniformity**






Image intensity distribution



Mixture of Gaussians model



Grey matter

White matter

Cerebrospinal fluid

## Segmentation into Different Tissues

Grey matter



White matter




Cerebrospinal fluid



Template tissue probability maps





**Normalization**



# [Preprocessing of sMRI]



- SPM25 [\[https://github.com/spm/spm\]](https://github.com/spm/spm)





## GUI of the SPM Toolbox


# Input




**T1-weighted and T2-weighted sMRI**

# Output



Correction for intensity non-uniformity




Intensity non-uniformity



**Correction of the T1-weighted Image for Intensity Non-uniformity**




Intensity non-uniformity




**Correction of the T2-weighted Image for Intensity Non-uniformity**

# Output

Segmentation



Segmentation









Image intensity distribution




Mixture of Gaussians model






**Tissue Classification based on a Mixture of Gaussians**

# Output

Normalization



Normalization




Grey matter






White matter



Cerebrospinal fluid




Template tissue probability maps





**Unified Segmentation and Normalization**

# Confirmation

Individual's normalized brain





MNI152 template brain



# Voxel-based Morphometry



- Without defining boundaries and modeling cortical surfaces
- CAT12 [<https://github.com/ChristianGaser/cat12>]






## CAT12 Toolbox

- Grey matter volume
  - Computed by multiplying voxel-wise grey matter probability by voxel volume
  - For a grey matter probability map in the native space or its modulated one in the standard space




**Voxel size:**  $1.5 \text{ mm} \times 1.5 \text{ mm} \times 1.5 \text{ mm}$   
**Voxel volume:**  $3.375 \text{ mm}^3$


## Computation of Grey Matter Volume for a Voxel or an Area



## Normalization and Modulation




## Features of Grey Matter Volume




AAL atlas




Hammers atlas




Harvard-Oxford atlas



Yeo atlas




AICHA atlas




Craddock atlas


## Brain Atlases



246 areas  
Brainnetome atlas






333 areas  
Resting-State Correlations atlas



360 areas  
HCP MMP 1.0 atlas

## Higher-resolution Brain Atlases

Atlas	# of regions	Horizontal	Sagittal	Coronal	Atlas	# of regions	Horizontal	Sagittal	Coronal
Hemispheric	2				Desikan	70			
Tissue	3				DKT	83			
Yeo-7	7				AAL	116			
Yeo-7-Lib	7				Glasser	360			
Yeo-17	17				CPAC200	200			
Yeo-17-Lib	17				Schaefer200	200			
HOS	21				Schaefer300	300			
Brodmann	41				Schaefer400	400			
HOC	48				Slab907	907			
JHU	48				Schaefer1000	1000			
PrincetonVis	49				Slab1068	1068			
PP264	58				Talairach	1105			


[<https://github.com/neurodata/neuroparc>]

## Repository of Human Brain Parcellations

# [Statistical Analysis of sMRI]

- Grey matter volume ~
  - Age +
  - Sex +
  - Education year +
  - Total intracranial volume (TIV) +
  - Memory performance

# Design matrix




# Output

Regression



Positive correlation




Negative correlation

- Input to machine learning models
  - Table of voxel-wise or area-wise grey matter volume values

	Features				
	Voxel or Area 1 grey matter volume	Voxel or Area 2 grey matter volume	Voxel or Area 3 grey matter volume	...	
Subject 1	-	-	-	-	-
Subject 2	-	-	-	-	-
Subject 3	-	-	-	-	-
:	-	-	-	-	-

- Grey matter volume map



[Adapted from Singh et al., 2020]

## Application of Deep Learning to Grey Matter Volume Maps

# Surface-based Morphometry

- Independent of registration and modulation
- Not applicable to subcortical regions
- FreeSurfer [\[https://github.com/freesurfer/freesurfer\]](https://github.com/freesurfer/freesurfer)




FreeSurfer


- Surface reconstruction
  - White matter surface: inner cortical boundary between the grey matter and white matter
  - Pial surface: outer cortical boundary between the grey matter and pia mater



White matter surface





Pial surface



[<https://www.physio-pedia.com/Meninges>]

## Cortical Surfaces beneath Cranial Meninges




Vertices and faces




[<https://surfer.nmr.mgh.harvard.edu/>]


## Surface Representation of the Brain



5124 vertices





8196 vertices




20484 vertices

## Surface Representation with Different Numbers of Vertices

- Cortical thickness
  - Distance between the inner (white matter surface) and outer (pial surface) cortical boundaries



[Gale and Huff, 2017]




## Features of Cortical Thickness

- Input to machine learning models
  - Table of vertex-wise or area-wise cortical thickness values

	Features				
	Vertex or Area 1 cortical thickness	Vertex or Area 2 cortical thickness	Vertex or Area 3 cortical thickness	...	
Subject 1	-	-	-	-	
Subject 2	-	-	-	-	
Subject 3	-	-	-	-	
:	-	-	-	-	

- Cortical thickness planar map




[Gao et al., 2022]

## Application of Deep Learning to Cortical Morphometric Maps

# Automated Brain Morphometry

- Employs deep learning algorithms to provide objective measurements of brain structures and their relationships
- Cortical metrics: Automated measurement of cortical thickness, surface area, and volume
- Subcortical segmentation: Precise delineation and volumetric analysis of deep brain structures



[Henschel et al., 2020; <https://github.com/Deep-MI/FastSurfer>]

## FastSurfer: Volumetric and Surface-based Cortical Thickness Analysis

Software	FDA Approved	MFDS Approved	EU CE Marked	URL	Country	Company	Character
NeuroQuant®	FDA approved	NA	CE marked	<a href="https://cortechs.ai/">https://cortechs.ai/</a>	USA	Cortechs Lab	NeuroQuant, LesionQuant*
Neuroreader®	FDA approved	NA	CE marked	<a href="https://brainreader.net/">https://brainreader.net/</a>	USA	Brainreader	
Icobrain	FDA approved	NA	CE marked	<a href="https://icomatrix.com/">https://icomatrix.com/</a>	EU	Icomatrix	Icobrain MS, Icobrain DM, Icobrain TBI, Icobrain ep
Quantib® Brain	FDA approved	NA	CE marked	<a href="https://www.quantib.com/">https://www.quantib.com/</a>	Netherlands	Quantib	
syMRI® Neuro	FDA approved	NA	CE marked	<a href="https://syntheticmr.com/">https://syntheticmr.com/</a>	Sweden		Multicontrast <sup>†</sup> Volumetry (GM, WM, CSF volume)
InBRAIN®	NA	MFDS approved	NA	<a href="https://www.inbrain.co.kr/">https://www.inbrain.co.kr/</a>	South Korea	Midas IT	
Neurol	NA	MFDS approved	NA	<a href="http://www.infomeditech.com/">http://www.infomeditech.com/</a> <a href="https://neurozen.ai/">https://neurozen.ai/</a>	South Korea	Infomeditech	
mdbrain	NA	NA	CE marked	<a href="https://www.qmenta.com">https://www.qmenta.com</a>	EU	Qmenta Inc	
DeepBrain®	NA	MFDS approved	CE marked	<a href="https://www.vuno.co/">https://www.vuno.co/</a>	South Korea	Vuno	
Atroscan	NA	MFDS approved	NA	<a href="http://jlkgroup.com">http://jlkgroup.com</a>	South Korea	JLK	

\*LesionQuant: measurement of WM lesion volume.

<sup>†</sup>Multicontrast: acquisition of multiple pre-defined contrast weighted images such as T1W, T2W and T2W FLAIR, as well as double inversion recovery and phase sensitive inversion recovery using synthetic MRI.

CSF = cerebrospinal fluid, EU CE Marked = conformite Europeenne marked, FDA = Food and Drug Administration, FLAIR = fluid attenuated inversion recovery, GM = grey matter, MFDS = Ministry of Food and Drug Safety, NA = not applicable, T1W = T1-weighted, T2W = T2-weighted, WM = white matter


[Lee et al., 2021]

## Commercially Available Brain Morphometry Software



[<https://www.cortechs.ai/>]



**NeuroQuant, Developed by Cortechs.ai, Offering Automated Volumetric Brain Structure Analysis**



neurophet AQUA



VUNO<sup>®</sup> Med-DeepBrain




NEUROZEN Neuro I


# sMRI as Anatomical Reference

- Anatomical alignment of other modalities of MRI
  - Within-subject between-modality registration

Rigid registration  
(global shift and rotation)




Registration




Within-subject within-modality

Affine registration  
(global shift, rotation, scale, and shear)



Registration



Within-subject between-modality

Deformable registration  
(local transformations)




Registration




Between-subject

[<https://kr.mathworks.com/help/medical-imaging/ug/medical-image-registration.html>]


## Image Registration



## Co-registration between sMRI and Functional MRI



Brain activity map overlaid on a functional image



Brain activity map overlaid on a structural image

## Anatomical Localization of Brain Activity